

CHANGE THE WAY YOU FEEL ABOUT TESTING

Introducing the **PELIKAN SUN™** electronic lancing device

TAKE THE PAIN OUT OF TESTING

The Pelikan Sun's patented electronic technology provides precision lancing control that virtually eliminates pain – a primary reason people don't test as often as they should – and significantly reduces finger trauma.

The Pelikan Sun's revolutionary lancing smoothness and precision means

- significant reduction of pain during lancing
- bleeding that stops almost immediately
- quick healing of the lanced site
- minimal residual bruising and sensitivity
- fingers that heal over time

NO TO PAIN, YES TO CONVENIENCE AND SAFETY


The Pelikan Sun uses a patented 50-lancet disk, ensuring users get a sterile, sharp lancet every time. Just pop in a new disk and you're ready to go. Never see or handle lancets again!



INTRODUCING FINGER-FRIENDLY™ TESTING

The Pelikan Sun offers 30 different levels of depth control, ensuring high sampling reliability while virtually eliminating pain.

Benefits of Pelikan Sun's patented technology:




Pelikan Sun

mechanical lancing

TAKE THE PAIN OUT OF TESTING

- Only Pelikan Sun instantaneously locates finger position and lances to an exact pre-set depth to ensure that the lancet reaches the capillary loops without penetrating nerve endings
- Pelikan Sun brings the lancet to a soft and controlled stop, eliminating the shock waves of pain caused by mechanical lancing devices



Pelikan Sun

mechanical lancing

SPONTANEITY OF BLOOD FLOW

- The slow and smooth retraction of the lancet from the skin allows blood to flow spontaneously, reducing the need to "milk" the site for blood



Pelikan Sun

mechanical lancing

CHANNEL QUALITY

- The precisely guided trajectory of Pelikan Sun's lancet ensures minimal trauma to the sample channel created by lancing
- Minimal residual bruising and sensitivity improves finger health over time
- Precise, electronically controlled lancing promotes immediate closure of the sample channel, so bleeding stops almost immediately

Taking the pain out of testing.

- USA: www.pelikantechnologies.com
1-877-PELIKAN
- Europe: www.pelikantechnologies.de
0800 735 7357
- Australia: www.diabetesnsw.com.au
1800 451 737



PELIKANSUN™